

GENDRON, INC.  
 520 W Mulberry St  
 Suite 100  
 Bryan Ohio 43506

CUSTOMER SERVICE  
 Phone: 800 / 537 - 2521  
 Fax: 419-636-9261  
 www.gendroninc.com

**Regency XL Wheelchairs Series 6500Q**  
**Seat Widths: 20" 22" 24" Weight Limit: 600 Pounds Price List Effective January 1 2009**

Date ..... Dealer PO #: .....  
 Dealer Name: ..... Contact: .....  
 Account Number: ..... Ship To: .....  
 Address: ..... Address: .....  
 City/State/Zip: ..... City/State/Zip: .....  
 Phone No.: (.....)..... Phone No.: (.....).....  
 Authorized Signature ..... Referral PO #: .....

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, fax or call-in order.

- | BASE WHEELCHAIR MODEL NUMBER  | OPTIONAL ACCESSORIES  |
|---|---|
| <input type="checkbox"/> 6500Q2.....\$1255  | <input type="checkbox"/> 12 oxygen holder..... \$135                                      |
| <b>Arm Style</b>  | <input type="checkbox"/> 44 8" semi pneumatic wheels.....standard                         |
| <input type="checkbox"/> permanent full length..... standard  | <input type="checkbox"/> 45 24" metal mag wheels solid tire..... \$350 pr                 |
| <b>Seat Dimensions W x D inches</b>   | <input type="checkbox"/> 59 leg guard..... \$100  |
| <input type="checkbox"/> 2 20 x 16 ..... standard   | <input type="checkbox"/> 75 solid padded back insert ..... \$160                          |
| <input type="checkbox"/> 3 22 x 16 ..... no charge  | <input type="checkbox"/> 80 safety belt auto buckle or velcro..... \$80                   |
| <input type="checkbox"/> 4 24 X 16..... no charge   | <input type="checkbox"/> 83 footplate heel loops ..... \$60                               |
| <b>Seat to Floor Distance</b>   | <input type="checkbox"/> 84 anti fold device ..... \$32                                   |
| <input type="checkbox"/> 20.5 inches..... standard  | <input type="checkbox"/> 85 wheeled anti tip attachment, front or rear ..... \$125        |
| <b>Front Riggings</b>   | <input type="checkbox"/> 86 wheel lock extensions ..... \$47                              |
| <input type="checkbox"/> 33 swing away foot rests ..... \$120                                       | <input type="checkbox"/> 93 patient records holder ..... \$55                             |
| <input type="checkbox"/> 35 swing away elevating leg rests..... \$235                               | <input type="checkbox"/> 94 solid padded seat insert ..... \$160                          |
| aluminum foot plates are standard   | <input type="checkbox"/> 95 IV rod / holder ..... \$105                                   |
| <b>Upholstery Colors</b>  | <input type="checkbox"/> 96 seat cushion, 2" foam..... \$120                              |
| <input type="checkbox"/> black <input type="checkbox"/> navy blue <input type="checkbox"/> burgundy | <input type="checkbox"/> 1294 oxygen tank/IV holder combination for pump mount..... \$240 |
| <input type="checkbox"/> gray <input type="checkbox"/> brown  | <input type="checkbox"/> 1295 oxygen tank/IV holder combination..... \$170                |
|   | <input type="checkbox"/> 1297 anti theft hoop/IV holder 20-24" wide..... \$148            |
|   | <input type="checkbox"/> 1299 anti theft hoop/IV/oxy. holder 20-24..... \$235             |
|   | <input type="checkbox"/> STNCL custom name stenciling..... \$75                           |
|   | Silk screen available. Contact customer service for pricing.                              |
|   | <b>Order Total \$ _____</b>   |

- Standard features include:**
- heavy duty textilene seat liner
  - 5/8" rear axle
  - rear axle and caster journal gussets
  - 1 1/4" square cross brace
  - seat strut supports
  - screw attached clothing guards
  - cushioned arm rests
  - through bolted seat upholstery
  - welded & reinforced caster fork stems
  - uni mold mag composite wheels

All specifications and dimensions are approximate.  
 Specifications subject to change without notice.  
 Gendron Inc. Form No. PLWCRegencyXL(0109)