

GENDRON, INC.
520 W Mulberry St
Suite 100
Bryan Ohio 43506

CUSTOMER SERVICE
Phone: 800 - 537 - 2521
Fax: 419-636-9261
www.gendroninc.com

Extra Care Alternating Pressure & Low Air Loss Support Surfaces
Weight Capacity: 350 to 1000 Pounds

Price List Effective January 1 2009

Date	Dealer PO #:
Dealer Name:	Contact:
Account Number:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone No.: (.....).....	Phone No.: (.....).....
Authorized Signature	Referral PO #:

To order, check boxes on the price list below. Add up prices, insert "total price", and mail, fax or call-in order.

EXTRA CARE SUPPORT SURFACES

- AL3680 36" x 80" alt. pressure & low air loss 350 lb. cap.....\$4000

EXTRA CARE BARIATRIC SUPPORT SURFACES

- AL3980 39" x 80" alt. pressure & low air loss 1000 lb. cap. \$5185
- AL4880 48" x 80" alt. pressure & low air loss 1000 lb. cap...\$5900
- AL5480 54" x 80" alt. pressure & low air loss 1000 lb. cap...\$6250
- ALD3980 39"x 80" digital AP/LAL 1000 lb. capacity\$5750
- ALD4880 48"x 80 " digital AL/LAL 1000 lb. capacity \$7250
- ALD5480 54"x 80" digital AL/LAL 1000 lb. capacity\$8400
- AP4880 48"x 80" alternating pressure 1000 lb. capacity.....\$4125
- AP5480 54"x 80" alternating pressure 1000 lb. capacity.....\$4450
- L4880 48"x 80" low air loss 1000 lb. capacity\$5175
- L5480 54"x 80" low air loss 1000 lb. capacity\$5550
- RL4880 rotation low air loss Z cell 48"x80" 1000 lb. cap.....\$8000
- RL5480 rotation low air loss Z cell 54"x 80" 1000 lb. cap.....\$8500

SPECIAL FEATURES

- All support surfaces and blowers made in USA
- All air cushions are urethane coated nylon
- All top sheets are low shear / low friction
- All top sheets are coated nylon and breathable
- CPR deflate utilizes quick disconnect connections
- Warranty: one year parts & labor for surface and electronics

All specifications and dimensions are approximate.
Specifications subject to change without notice.
Gendron Inc. Form PLSTSUPSUR (0109)