

GENDRON, INC
 520 W Mulberry St
 Suite 100
 Bryan Ohio 43506

CUSTOMER SERVICE
 Phone: 800-537-2521
 Fax: 419-636-9261
 www.gendroninc.com

Regency XLC Wheelchairs Custom Bariatric Manual and Power Wheelchairs
Seat Widths: 20" plus Weight Limit: 700 Pounds Plus Price List Effective January 1 2009

Date P.O. #:
 Dealer Name: Contact:
 Account Number: Ship To:
 Address: Address:
 City/State/Zip: City/State/Zip:
 Phone No.: (.....) Phone No.: (.....)
 Authorized Signature Therapist / Referral:

Please Note: all XLC wheelchairs are custom built and are non-returnable.

RIGID FRAME WHEELCHAIR BASE

- 7850 Manual base \$3500
- 7800 Power base \$9000
- Power base right control Power base left control.....N/C

Seat Width Please call for seat widths greater than 40"

- | | |
|---|---|
| <input type="checkbox"/> 20" standard | <input type="checkbox"/> 31" \$2100 |
| <input type="checkbox"/> 21" \$375 | <input type="checkbox"/> 32" \$2300 |
| <input type="checkbox"/> 22" \$500 | <input type="checkbox"/> 33" \$2475 |
| <input type="checkbox"/> 23" \$625 | <input type="checkbox"/> 34" \$2675 |
| <input type="checkbox"/> 24" \$750 | <input type="checkbox"/> 35" \$2850 |
| <input type="checkbox"/> 25" \$950 | <input type="checkbox"/> 36" \$3025 |
| <input type="checkbox"/> 26" \$1125 | <input type="checkbox"/> 37" \$3225 |
| <input type="checkbox"/> 27" \$1325 | <input type="checkbox"/> 38" \$3400 |
| <input type="checkbox"/> 28" \$1500 | <input type="checkbox"/> 39" \$3600 |
| <input type="checkbox"/> 29" \$1675 | <input type="checkbox"/> 40" \$3775 |
| <input type="checkbox"/> 30" \$1925 | call for adjustable seat width option |

Note: during manufacture rear wheels on power base chairs can be "tucked" under seat frame. Consult factory for narrowest overall width based on seat dimension.

Seat Construction

- solid pan, no cushion, metal deck \$t
- solid with 2" foam \$225
- solid with 4" foam \$275
- Profiler XL bariatric pressure relief support cushion call

Back Style

- permanent standard
- manual recliner 45 degrees. 60 degrees \$2500
- removable headrest extension 12" high \$375

Power Accessories / Back / Tilt in Space

- power recliner 45 degrees. 60 degrees \$3500
- power recliner 90 degrees \$4000
- power tilt in space 30 degrees. 45 degrees \$4900

Back Support Surface

- solid back, no padding standard
- solid, 2" foam pad \$225
- solid, 4" foam pad \$275

Arm Style

- full length, height adjustable, 7"-9.5", removable standard
- full length, height adjustable, 9" - 11", removable \$275

Front Riggings/Heavy Duty/Reinforced

- footrest for seat height 14"-16" \$625
- footrest for seat height 16"-20" \$500
- elev. legrests for seat height 14"-16" \$675
- elev. legrests for seat height 16"-20" \$575
- power articulating elevating leg rests \$3125
- extra wide calf pad supports, state width \$200 pr

Frame Colors

- black standard

Seat depth

- | | |
|---|---|
| <input type="checkbox"/> 18" standard | <input type="checkbox"/> 23" \$800 |
| <input type="checkbox"/> 19" \$575 | <input type="checkbox"/> 24" \$1000 |
| <input type="checkbox"/> 20" \$575 | <input type="checkbox"/> 25" \$1175 |
| <input type="checkbox"/> 21" \$575 | <input type="checkbox"/> 26" \$1375 |
| <input type="checkbox"/> 22" \$800 | |

Seat Height*

- 14" 15" avail. on manual wheelchairs only \$950
- 16", 17", 18", 19" manual or power wheelchairs \$675
- 20" manual or power wheelchairs \$325
- 21" manual or power wheelchairs standard

*Measured from floor to seat pan without foam padding.

For any other seat heights/depths not listed, please call customer service

Seat Plane / Angle

- 3 degrees standard
- other seat angles, specify distance from floor for:
 front edge of seat _____ rear edge of seat \$375

Back Height

- 16" standard
- 17" 18" 19" 20" \$250
- 21" 22" 23" 24" \$375
- 25" 26" 27" 28" \$500

Back Angle

- back angle 7 degrees standard
- other back angle, specify _____ \$375

OPTIONAL ACCESSORIES

- 12 oxygen tank holder. (see Note) \$200
- 80 safety belt/auto buckle or velcro closure \$75
- 83 heel loops \$50
- 85 wheeled anti-tip brackets / rear standard
- 86 wheel lock extensions, manual chairs \$50
- 95 IV rod and holder (see Note) \$200
- 1295 oxygen tank/IV rod holder combination (see Note) \$350
- P20529 battery 12V, gell cell \$200 ea

Note: IV / oxygen holders not available on recliners & tilt in space

Order Total \$ _____

All specifications and dimensions are approximate. Specifications subject to change without notice. Gendron Inc. Form No. PLWC RegencyXLC(0109)

