

GENDRON, INC.  
 520 W Mulberry St  
 Suite 100  
 Bryan Ohio 43506

CUSTOMER SERVICE  
 Phone: 800.537.2521  
 Fax: 419-636-9261  
 www.gendroninc.com

**EC + Extra Care Bariatric Patient Lift  
 Weight Limit: 700 – 1000 Pounds**

**Price List Effective January 1 2009**

Date .....  
 Dealer Name: .....  
 Account Number: .....  
 Address: .....  
 City/State/Zip: .....  
 Phone No.: (.....) .....  
 Authorized Signature .....

P.O. #: .....  
 Contact: .....  
 Ship To: .....  
 Address: .....  
 City/State/Zip: .....  
 Phone No.: (.....) .....  
 Therapist / Referral: .....

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, [e-mail](#), fax or call-in order.

**BASE MODEL NUMBER**

- GT700 700 pound capacity ..... \$7900
- GT700S 700 pound capacity with digital scale ..... \$9150
- GT1000 1000 pound capacity ..... \$13800
- GT1000S 1000 pound capacity with digital scale ..... \$14950

**Control / Power Source**

- Battery operated, pendant hand control ..... standard
- Battery operated, integral thumb push button..... standard

**Casters**

- Four total, two rear locking ..... standard

**Lift Function & Suspension**

- Four point lift padded bar ..... standard
- Fully electric lift operation ..... standard

**Accessories & Slings**

- GT3XL bariatric sling (400 – 700 lbs.)..... standard
- GT4XL bariatric sling (400 – 800 lbs.)..... \$385
- GT5XL bariatric sling (400 – 1000 lbs.)..... \$480
- GT10151-1 stretcher frame & sling..... \$1400
- GT700B battery for models GT700 / GT700S..... \$450
- GT1000B battery for models GT1000 / GT1000S ..... \$450
- GT700CH battery charger for models GT700 / GT700S ..... \$570
- GT1000CH battery charger for models GT1000 / GT1000S. . \$570

**STANDARD FEATURES INCLUDE:**

- fully electric lift operation
- two easy install batteries
- battery charger
- size 3XL bariatric sling (400 – 700 lbs.)
- pendant control and integral thumb push button controls
- manually retractable base legs for stability
- four point sling suspension
- lift range: 31"-74" (model GT700) / 39" – 76" (model GT1000)
- base size opened: model GT700 / 44.5" model GT1000 / 60"
- base size closed: model GT700 / 30" model GT1000 / 36"
- rear locking casters
- under bed clearance: model GT700 / 4.5" model GT1000 / 6.75"
- Warranty: five years welded structure / one year electrical

**Order Total \$** \_\_\_\_\_

All specifications and dimensions are approximate  
 Specifications subject to change without notice.  
 Gendron Inc. Form No. BARLIFT (0109)