

GENDRON, INC.
 520 W Mulberry St
 Suite 100
 Bryan Ohio 43506

CUSTOMER SERVICE
 Phone: 800.537.2521
 Fax: 419-636-9261
 www.gendroninc.com

Price List Effective January 1 2009

Regency XL 2000 Bariatric Manual Wheelchairs
Seat Widths: 26" 28" 30" 32" 34"

Model Series: 6800
Weight Limit: 850 Pounds

Date
 Dealer Name:
 Account Number:
 Address:
 City/State/Zip:
 Phone No.: (.....)
 Authorized Signature

P.O. #:
 Contact:
 Ship To:
 Address:
 City/State/Zip:
 Phone No.: (.....)
 Therapist / Referral:

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, [e-mail](#), fax or call-in order.

BASE WHEELCHAIR MODEL NUMBER

68 base wheelchair \$3265

Seat Dimensions W x D inches

- 6826-20 26 X 20 standard
- 6828-20 28 X 20 \$200
- 6830-20 30 X 20 \$325
- 6832-20 32 x 20 \$450
- 6834-20 34 x 20 \$550
- 6826-22 26 X 22 \$250
- 6828-22 28 X 22 \$375
- 6830-22 30 X 22 \$500
- 6832-22 32 x 22 \$625
- 6834-22 34 x 22 \$725

Seat to Floor Height

- (8) 18" (1) adult N/C
- (7) 17" (1) adult N/C
- (6) 16" (1) hemi N/C
- (5) 15" (2) super hemi (includes 20" and 6" wheels) \$320
 - (1) 24" rear / 8" caster standard for 16" 17" 18" seat hts.
 - (2) 20" rear / 6" caster required for 15" seat ht.

Arm Style

- (-62) full length, removable standard
- (-63) full length, removable 8" - 12" height adjustable \$150

Push Bar

- Full width removable push bar standard

Back Style

- adjustable height 16", 17" & 18" (factory set at 16") standard
- (-78) power assist 45 degree reclining back, 36" high \$1100
 note: push bar and rear anti tips are standard on recliners

Frame Colors

- Black standard

Front Riggings

- (-26) adjusta front footrests \$240
- (-27) adjusta front elevating legrests \$420
- (-27A) adjusta front articulating elevating legrests \$525
- custom length extension tubes _____ Inches \$185

Wheel Options

- 24" metal wheels solid tires for 16" 17" & 18" height standard
- 20" mag wheels with semi pneumatic 15" ht. standard
 Note: wheels have steel hand rims

Caster Options

- 8" x 2" semi pneumatic wheels for 16" 17" 18" height standard
- 6" x 2" semi pneumatic wheels for 15" height standard

IV & Oxygen Tank Holder

- (-12) oxygen holder \$135
- (-13) oxygen holder for recliner \$135
- (-14) IV rod / holder for recliner \$105
- (-95) IV holder height adjustable \$105
- (-1295) oxygen / IV holder combination \$170

Miscellaneous Accessories

- (-59) leg guard \$100
- (-80) safety belt, specify auto buckle Velcro. \$80
- (-85) anti tippers, rear with wheels \$125
- (-86) wheel lock extensions \$50

Standard features include:

- heavy duty steel frame, 5/8" rear axle
- double cross-braces
- cushioned arm rests with steel clothing guards
- through bolted seat upholstery
- black nylon upholstery with heavy duty seat liner
- caster fork & caster journal-gussets reinforced & welded
- back angle 10 degrees
- 8" x 2" semi-pneumatic caster wheels
- 24" metal mag wheels with steel hand rims
- adjustable seat to floor height
- all models feature heavy-duty adjusta front riggings
- full length push bar across back

All specifications and dimensions are approximate
 Specifications subject to change without notice
 Gendron Inc. Form
 No.: PLWC 6800 Series (0109)