

GENDRON, INC.
 400 E LUGBILL RD
 ARCHBOLD OHIO 43502

CUSTOMER SERVICE
 Phone: 800.537.2521
 Fax: 419.446.2631
 website: <http://www.gendroninc.com>
 e-mail: sales@gendroninc.com

**EC + Extra Care Bariatric Patient Lift
 Weight Limit: 700 – 1000 Pounds**

Price List Effective March 1 2006

Date.....
 Dealer Name:
 Account Number:
 Address:
 City/State/Zip:
 Phone No.: (.....).....
 Authorized Signature

P.O. #:
 Contact:
 Ship To:
 Address:
 City/State/Zip:
 Phone No.: (.....).....
 Therapist / Referral:

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, [e-mail](mailto:sales@gendroninc.com), fax or call-in order.

BASE MODEL NUMBER

- GT700 700 pound capacity..... \$8,160
- GT700S 700 pound capacity with digital scale \$9,590
- GT1000 1000 pound capacity \$14,050
- GT1000S 1000 pound capacity with digital scale \$15,500

Control / Power Source

- Battery operated, pendant hand control..... standard
- Battery operated, integral thumb push button..... standard

Casters

- Four total, two rear locking standard

Lift Function & Suspension

- Four point lift padded bar standard
- Fully electric lift operation standard

Accessories & Slings

- GT3XL bariatric sling (400 – 700 lbs) standard
- GT4XL bariatric sling (400 – 800 lbs.) \$385
- GT5XL bariatric sling (400 – 1000 lbs.) \$480
- GT10151-1 stretcher frame & sling..... \$1110
- GT700B battery for models GT700 / GT700S \$280
- GT1000B battery for models GT1000 / GT1000S \$360
- GT700CH battery charger for models GT700 / GT700S \$280
- GT1000CH battery charger for models GT1000 / GT1000S. . \$360

STANDARD FEATURES INCLUDE:

- fully electric lift operation
- two easy install batteries
- battery charger
- size 3XL bariatric sling (400 – 700 lbs.)
- pendant control and integral thumb push button controls
- manually retractable base legs for stability
- four point sling suspension
- lift range: 31"-74" (model GT700) / 39" – 76" (model GT1000)
- base size opened: model GT700 / 44.5" model GT1000 / 60"
- base size closed: model GT700 / 30" model GT1000 / 36"
- rear locking casters
- under bed clearance: model GT700 / 4.5" model GT1000 / 6.75"

Order Total \$ _____

All specifications and dimensions are approximate
 Specifications subject to change without notice.
 Gendron Inc. Form No. BARLIFT (0306)